

Centrifuge Tubes

Features

- High quality, ultra-transparent polypropylene
- Sterile, DNase/RNase free, non-pyrogenic
- Maximum centrifugal force to 9500×g RCF (50ml tubes) and 8000×g RCF (15ml tubes)
- Withstand temperatures of -80°C and -20°C
- Large white writing area
- Leak-proof seal design, single-hand operation

i-labware
more than just equipment



Ordering Information

Cat No.	Description	Package Quantity
C52611	15ml centrifuge tubes, PP, sterile	25/bag, 20bags/case
C52621	15ml centrifuge tubes, PP, sterile, with shelves	50/rack, 10 racks/case
C52811	50ml centrifuge tubes, PP, sterile	25/bag, 20 bags/case
C52821	50ml centrifuge tubes, PP, sterile, with shelves	25/rack, 20 racks/case