

i-Mix 3D Rotating Mixer

Cat. No.
E30010

i-Labware
more than just equipment



Equipment Specifications

Time Range	1s - 999 mins
Reciprocating Rotation Degree	20° - 99°
Power	AC100/230V, 50/60Hz, ≤4W
Fuse	0.5A /1.0A
Speed	10-40rpm
Dimension	240 x 145 x 210 mm

Features

- Can be operated in variety of rotation modes by horizontal, vertical or any angled rotation
- Compact size. Easily moved from benchtop to refrigerator or incubator

Includes 3 different rotation frames:

- 1. 5 x 10ml /15 ml (Φ15 mm-Φ16 mm) centrifuge tubes + 6 x 5ml /7 ml (Φ12 mm-Φ13 mm) centrifuge tubes
- 2. 3 x 50 ml centrifuge tubes + 18 x 0.5 ml centrifuge tubes
- 3. 18 x 1.5/2 ml centrifuge tubes