

i-Mix Smart Vortexer

Cat. No.
E30011

i-Labware
more than just equipment



Features

- Convenient and easy handling of touch-to-mix (short mix) mode
- Elegant and compact design
- Outer surfaces are made of chemical-resistant plastic with a thick base for steady operation
- Adjustable speed to accommodate to different experimental requirements
- Built-in, portable 12V low voltage power adapter

Equipment Specifications

Mixing method	Circular
Shaking Orbit	4.5 mm
Speed	0 - 2800 rpm
Power	AC90-230V, 15W
Dimensions	110 × 110 × 85 mm
Net Weight	1.2 kg