

i-Vac Vacuum Pump

Cat. No.
E40010

i-Labware
more than just equipment



Features

- Oil-free working cavity
- Corrosion-resistant instrument parts (polyphenylene sulfide & aluminium alloy)
- Pressure indicator and negative pressure gauge for convenient operation
- Small footprint. Easy to move about

Equipment Specifications

Rate of Flow	26L/min
Vacuum Degree	0.08 Mpa
Negative Pressure Adjust Range	0-0.08 Mpa
Filtration Bottle Capacity	250 ml
Power	AC 220V/110V, 50HZ/60HZ, 60W