

PCR Plates & Accessories

Features

- USP class VI certified polypropylene plates
- Thin-walled for optimal heat transfer
- Free of detectable DNase, RNase, DNA, PCR Inhibitors, ATP, Endotoxins
- Designed to minimize sample loss
- Broad compatibility with RT-PCR instruments

i-labware
more than just equipment



Ordering Information

Cat No.	Description	Package Quantity
C20010	0.2ml PCR 96 well plate, clear, non-skirted	10/bag, 5 bags/case
C20020	0.1ml qPCR 96 well plate, clear, sub-semi skirted	10/bag, 5 bags/case
C20021	0.1ml qPCR 96 well plate, low profile, white, semi skirted	10/bag, 5 bags/case
C20022	0.2ml qPCR 96 well plate, clear, sub-semi skirted	10/bag, 5 bags/case
C20023	0.1ml qPCR 96 well plate, low profile, clear, semi skirted	10/bag, 5 bags/case
C30010	Storage sealing membrane	10/bag, 10 bags/case
C30011	Aluminum sealing membrane	10/bag, 10 bags/case
C30020	PCR 96 well plate sealing membrane	10/bag, 10 bags/case
C30021	qPCR 96 well plate optical sealing membrane	10/bag, 10 bags/case

i-Labware Pte Ltd

enquiry@i-labware.com | www.i-labware.com

Instrument Compatibility for PCR Plate

	Model	C20010	C20020	C20021	C20022	C20023
	Standard Well (0.2ml, height 21mm)	✓			✓	
	Low Profile Well (0.1ml, height 16mm)		✓	✓		✓
Applied Biosystems	7000				✓	
	7300				✓	
	7500				✓	
	7900HT Fast (with tray) (0.1ml type)		✓			
	7900HT (with tray) (0.2ml type)				✓	
	StepOne Plus		✓			
	7500Fast		✓			
	ViiA7 Fast (0.1ml model)		✓			
	ViiA7 (0.2ml model)				✓	
	QuantStudio Fast 96-well Block (0.1ml model)		✓			
	QuantStudio 96-well Block (0.2ml model)				✓	
	5700				✓	
	7700				✓	
	310/3100/3130				✓	
	3700/3730				✓	
Roche	LightCycler 480			✓		
	LightCycler 480 II			✓		
	LightCycler 96			✓		
Bio-Rad	iCycler iQ	✓				
	iQ5	✓				
	MyiQ	✓				
	Chromo4	✓				
	Opticon I			✓		✓
	Opticon II			✓		✓
	CFX96			✓		✓
Others	Stratagene Mx4000/3000P/3005P	✓				
	Eppendorf MasterCycle ep Realplex	✓				
	Biometra Topical Thermocycler	✓				