

i-Cube Microplate Incubator

Cat. No.
E10030

i-Labware
more than just equipment



Features

- Microprocessor ensures accurate time and temperature control
- Alarm signals upon operation completion
- User-friendly interface displays setting and actual temperatures simultaneously



Equipment Specifications

Temperature Range	R.T.+5°C - 70°C	Heating Time	≤25 min (R.T. to 70°C)
Timing	1min - 99h59min	Capacity	2 Plates
Temperature Precision	≤ ± 0.5°C	Power	150W
Display Accuracy	0.1°C	Dimension	248 x 264 x 157 mm
Module Temperature Uniformity	≤ ± 0.5°C	Net Weight	3.0 kg

i-Labware Pte Ltd

enquiry@i-labware.com | www.i-labware.com