

i-Spin Mini Centrifuge

Cat. No.
E20010/
E20011

Features

- Fits 8 tubes compared to the usual 6, allowing for 1.5x greater centrifugal force at the same speed
- Convenient interchangeable rotors with easy clip-on fix
- Includes rotor accessories and tube adaptors.
- Equipped with 2 rotor types
 1. 8 x 1.5ml rotor (suitable for 8 x 1.5ml / 2.0ml tubes with additional 0.2ml & 0.5ml tube adapters)
 2. 8 x 0.2ml x 2-strip rotor (suitable for 8 x 0.2ml PCR tubes)

i-Labware
more than just equipment



Equipment Specifications	i-Spin 6k Mini Centrifuge	i-Spin 10k+ Mini Centrifuge
Cat No.	E20010	E20011
Speed	6000 rpm	3000 - 10 000 rpm
Centrifugal Force	3000 g	1000 - 7500 g
Time Range	1s - 99min59s	1s - 999 mins
Noise	≤45 dB	≤55 dB
Power	AC110/220V, 50/60Hz, ≤45W	AC110/220V, 50/60Hz, ≤45W
Fuse	0.5A /1.0A	0.5A /1.0A
Dimension	176 x 156 x 121 mm	176 x 156 x 121 mm
Net Weight	1.5 kg	1.5 kg